

PROJECT WAS DISCUSSED WITH MCTP-STAFF

metabolomics-projects@uni-heidelberg.de

ONLINE REQUEST SUBMITTED (PER INDIVIDUAL ANALYSIS)

hmls.corefacilities.org/

SAMPLES PREPARED ACCORDING TO USER GUIDELINES

- Initials + consecutive numbering as label, e.g. MM_01
 - No adhesives only written labels
 - Sample sets >15 samples, please randomize sample order
- For correct sample weights/volumes, please check user guidelines
 - Solids: Weighed and pulverized in 1.5/2ml tubes
 - Solids whole: weighed sample in 2ml tubes
 - Liquids: Not overfilled in 1.5/2ml tubes

SAMPLE HANDOVER

- Transfer on dry ice
- Delivery in plastic or paper box

Submission times:

MON 10:30 – 11:30, TUE 9:00 – 11:00, WED 13:30 – 15:00, THU 9:00 – 11:00; other times upon prior consultation via metabolomics-samples@uni-heidelberg.de

Platform